

Z327 Universal Series



Tissue Culture Pkg.

Blood Tube Pkg.

As low as \$8,600

- 2 bundle package choices:

Tissue Culture: up to 8 x 50ml and 16 x 15ml conical tubes

Blood Tube: up to 28 x 8-10ml and 32 x 1.8-7ml blood tubes

- Extremely versatile, 17 rotor options
- Refrigerated model available (Z327-K)
- High speed rotors up to **23,542 xg**

BUNDLE ORDERING INFORMATION*

		Reg.	Bundle
Z327-BND	Z327 Ambient Tissue Culture Bundle	\$9,867	\$8,600
Z327-BND-BD	Z327 Ambient Blood Tube Bundle	\$10,096	\$8,600
Z327-K-BND	Z327-K Refrigerated Tissue Culture Bundle	\$12,952	\$11,200
Z327-K-BND-BD	Z327-K Refrigerated Blood Tube Bundle	\$13,181	\$11,200

Z366 Universal Series



Z366-K

As low as \$9,000

- Max. swing-out capacity **16 x 50ml** or **40 x 15ml**
- High speed rotors up to **24,000 xg**
- Refrigerated model available (Z366-K)
- Capacity up to 6 x 250ml

BUNDLE ORDERING INFORMATION*

		Reg.	Bundle
Z366-BND	Z366 Ambient Tissue Culture Bundle	\$11,178	\$9,000
Z366-K-BND	Z366-K Refrigerated Tissue Culture Bundle	\$14,949	\$11,600

Z446 Universal Series



Z446-K

As low as \$10,700

Offers valid until June 30th 2024

- Max. swing-out capacity:
28 x 50ml, 68 x 15ml, and 20 x microplates
- Temp Range -20°C to +40°C (Z446-K)
- Bundle options for **tubes** or **plates!**
- Truly universal, **14** rotor options



BUNDLE ORDERING INFORMATION*

		Reg.	Bundle
Z446-BND	Z446 Ambient Tissue Culture Bundle	\$12,980	\$11,200
Z446-MPBND	Z446 Ambient Microplate Bundle	\$13,151	\$10,700
Z446-K-BND	Z446-K Refrigerated Tissue Culture Bundle	\$16,885	\$14,500
Z446-K-MPBND	Z446-K Refrigerated Microplate Bundle	\$17,055	\$14,100

2024 Hermle Centrifuge Bundles



Z206 Compact Centrifuge

- Includes 6 x 50ml rotor
- Compact, only 11 inches wide
- Optional rotors, for tubes from 0.2ml to 50ml

BUNDLE ORDERING INFORMATION

Z206-BND Z206-A Compact centrifuge with 6 x 50ml rotor ~~Reg. \$3,670~~ **Bundle \$2,800**



Z307 Deepwell Microplate

- Swing out plate rotor accepts deepwell plates or can house up to 3 standard plates on each arm
- High Speed microtube rotors (up to 18,624xg)
- Fixed angle rotors for 0.2ml to 50ml

BUNDLE ORDERING INFORMATION*

Z307-BND-MP Z307 with Microplate rotor ~~Reg. \$7,246~~ **Bundle \$5,700**



Z216-MK Refrigerated Microcentrifuge

- Temp Range -20°C to +40°C
- 9 rotor options for 0.2ml to 5ml tubes
- High speed, up to 21,631xg
- 2 innovative rotor bundle options!



BUNDLE ORDERING INFORMATION*

Z216-K-BND44 Z216-MK Refrigerated w/ 44 place rotor ~~Reg. \$7,825~~ **Bundle \$6,400**
 Z216-K-CMB Z216-MK Refrigerated w/ COMBI-Rotor™ ~~Reg. \$8,506~~ **Bundle \$6,900**

Ambient versions also available



Offers valid until June 30th 2024

Additional Bundles Available!*

* inquire for models not shown

